

3 in 1

Body「Contouring」 & Laser「Therapy」

Microwaves at 2450MHz



Shallow handpiece



Deep handpiece




685nm Laser Physical
Therapy Handle



Working Theory

FREEZEWAVE is a non-invasive microwave therapy that selectively targets fat cells using controlled microwave energy to safely, effectively, and non-invasively reduce localized fat deposits and cellulite. Microwave stimulation causes fat cells to release lipids, disrupts the connective tissue responsible for the orange peel effect, and simultaneously stimulates collagen and elastin production in surrounding tissues. This helps improve skin firmness and elasticity, resulting in smoother skin.

The device operates at a frequency of 2.45GHz, the energy passes directly through the surface layers without overheating them. Only 20% of the cold wave energy produces a thermal effect, while the majority (approximately 80%) effectively targets subcutaneous fat, disrupting fat cells and dissolving localized fat deposits.




685nm Laser Physical Therapy

Used for superficial areas


It combines bio stimulation and optomechanical stimulation to provide powerful and non-addictive pain management while promoting tissue healing. Based on the principle of laser thermal action, the laser with precisely controlled energy directly irradiates the skin and stimulates the deep muscle tissue, which can cause a series of biological effects, promote blood circulation, improve body function, and at the same time play a role in soothing and activating the muscles and anti-inflammatory, anti-swelling and pain relief. It can also help break down fat.




FREEZEWAVE




A




B




C




A



B



A



B

Histological images of tissue with human abdominal fat. (A): Control (magnification x10). (B): Sample from the same patient immediately after treatment with FREEZEWAVE (magnification x20). Image (B) clearly shows the ruptured lipocytes and initial hyperaemia with dilatation of the blood vessels.

80%

of Freezewave energy
reaches fat layer to
destroy fat cells

VS

20%

superficial heat
retained in the
dermis





Product Parameter

【Fuse】	10A
【Power】	800W
【Screen】	12.1-inch touch screen
【Voltage】	AC110V-220V
【Frequency】	2450MHz
【Wavelength】	685nm red Laser
【Output Power】	12J~82J
【Work handle】	3 handpieces (Shallow/Deep/Laser therapy)
【Packing size】	113cm*63cm*61 cm
【N. W./G. W. 】	32KGS / 54KGS



Parts Introduction



Two colors to choose



Grayish White



Grayish Black

12.1-inch touch screen



Smart UI With, Easy to Operate





Easy operation



Quick response






Shallow Handpiece

Treatment area 5.2*5.2cm

Applicable Areas: Skin Tightening, Double chin, chin fat, inner arms, superficial abdominal fat etc.

5.2cm


A white, ergonomic handpiece with a silver-colored circular treatment head. A green ring is visible on the handle. A red dimension line indicates the diameter of the treatment head is 5.2cm.

Deep Handpiece

Treatment area 6.5*6.5cm

Applicable Areas: Body slimming, Thick fat areas such as the waist, abdomen, thighs, buttocks, and sides, as well as stubborn fat areas.

6.5cm

A white, ergonomic handpiece with a silver-colored circular treatment head. A green ring is visible on the handle. A red dimension line indicates the diameter of the treatment head is 6.5cm.

Laser Therapy

685nm red light Therapy

Treatment area 6.7*6.7cm

Applicable Areas: conditions such as frozen shoulder, cervical spondylosis, lumbar muscle strain, acupuncture point irradiation and Shallow fat removal.

6.7cm

A white, ergonomic handpiece with a silver-colored circular treatment head. The treatment head features a grid of red dots. A red ring is visible on the handle. A red dimension line indicates the diameter of the treatment head is 6.7cm.

REEZEWAVE

Deep Reduced Fat

Microwaves at 2450MHz

Painless body shaping



Reduced Fat



Tighten Skin



Body shaping

Skin Tightening



Microwaves at 2450MHz

Improve skin sagging

Remove double chin



Laser Physical Therapy

Physical Therapy

Reduce fat,slimming

685nm red light Therapy

frozen shoulder,neck pain



685nm Laser Therapy



Reduced Fat & Skin Tightening & Laser Therapy

FREEZEWAVE



SHALLOW HANDPIECE

Treatment area 5.2*5.2cm



Tighten
skin



Improve
sagging



Eliminate
cellulite



DEEP HANDPIECE

Treatment area 6.5*6.5cm



Under skin
reach



Reduce
fat



Body
shaping



LASER THERAPY HANDPIECE

Treatment area 6.7*6.7cm

685nm red light Therapy



Physical
Therapy



lumbar
muscle strain



Reduce
fat



「2.45GHZ」 MICROWAVE FREEZEWAVE




Shallow handpiece



Deep handpiece



685nm Laser
Physical Therapy
Handle



Skin Tightening

Uses microwave energy to activate collagen and fat layers for better lifting and tightening.

Double Chin Fix

Reduces sagging double chin and defines your jawline.

Reduction Fat

Effectively reduces stubborn fat in areas like arms, belly, and legs.

Laser Therapy

For conditions such as frozen shoulder, cervical spondylosis, lumbar muscle strain, acupuncture point irradiation skin beauty and Helps break down fat.

3in1 Freezewave Function

Remove
abdomen fat



Remove
back fat



Remove
arm fat



Remove
thigh fat



Skin tightening



Double chin



Frozen
shoulder



lumbar strain





ADVANTAGE



Microwaves at 2450 Mhz
Frequency, 800W high-power
output energy



12.1-inch touch screen,
Intelligent system, easy to operate



80% energy acts on fat,
20% acts on surface heat effect



Function: dissolves fat, Skin Tightening,
Laser therapy: frozen shoulder,
lumbar strain, fat breakdown



Instrument Use Case

FREEZE WAVE





Before & After

FREEZE WAVE



Product Function



Product Advantage

- Microwaves at 2450 Mhz Frequency, 800W high-power output energy.
- 12.1-inch touch screen, Intelligent system, easy to operate
- Function: Intelligent Freezewave therapy dissolves fat, Localized Fat Deposits, Skin Tightening, 685nm laser therapy: frozen shoulder, lumbar strain, fat breakdown.
- 80% Freezewave energy targeting on adipocytes.
- 20% Secondary superficial thermal effect.

Before and After

