

4 in 1



Body Contouring & Laser Therapy Machine



**Pocket
handpiece**



**Shallow
handpiece**



**Deep
handpiece**




**685nm Laser
Physical Therapy**

Working Theory


FREEZEWAVE is a non-invasive microwave therapy that selectively targets fat cells using controlled microwave energy to safely, effectively, and non-invasively reduce localized fat deposits and cellulite. Microwave stimulation causes fat cells to release lipids, disrupts the connective tissue responsible for the orange peel effect, and simultaneously stimulates collagen and elastin production in surrounding tissues. This helps improve skin firmness and elasticity, resulting in smoother skin.

The device operates at a frequency of 2.45GHz, the energy passes directly through the surface layers without overheating them. Only 20% of the cold wave energy produces a thermal effect, while the majority (approximately 80%) effectively targets subcutaneous fat, disrupting fat cells and dissolving localized fat deposits.



685nm Laser Physical Therapy

It combines bio stimulation and optomechanical stimulation to provide powerful and non-addictive pain management while promoting tissue healing. Based on the principle of laser thermal action, the laser with precisely controlled energy directly irradiates the skin and stimulates the deep muscle tissue, which can cause a series of biological effects, promote blood circulation, improve body function, and at the same time play a role in soothing and activating the muscles and anti-inflammatory, anti-swelling and pain relief. It can also help break down fat.





Product Parameter

- Fuse:10A
- Power:1000W
- Screen:12.1 inch
- Frequency:2450MHz
- Laser Wavelength : 685nm
- Cooling Temperature : 2°C~10°C
- Water tank capacity : 6L
- Water tank volume : 216* 240*119mm
- Output Power : 12~82J (body); 6~40J (face)
- Voltage : AC110V-220V
- Packing size : 117cm*65cm*60cm
- N.W. /G.W. : 40kg / 65kg



Parts Introduction



12.1-inch touch screen



Smart UI With, Easy to Operate







Easy operation



Quick response





Shallow Handpiece

Adopting ergonomic design, intelligently enhancing the shaping of cold waves in various parts of the body , imultaneously protect the skin, ensuring safety and comfort.

2.45 GHz Ultra-Microwave Technology

- Can precisely target the subcutaneous fat layer at a focal depth of 3-12mm
- Rapidly penetrates the epidermis and dermis, delivering more concentrated thermal energy.

Treatment Start Button



Cooling temperature:
-10°C

Protects the epidermis and dermis, and controls thermal stimulation ●

4 HANDPIECES FOR COMPLETE BODY TREATMENT

Pocket Handpiece

Focal Depth 3mm

- Treatment area 4.2*4.2cm
- Applicable Areas: Full face lifting, neck, jawline, around the eyes, around the mouth, and other delicate areas.

4.2cm

A close-up of the Pocket Handpiece, showing a circular treatment area with a diameter of 4.2cm. The handpiece is grey with a green ring around the treatment area.

Shallow Handpiece

Focal Depth 7mm

- Treatment area 5.2*5.2cm
- Applicable Areas: Skin Tightening, Double chin, chin fat, inner arms, superficial abdominal fat etc.

5.2cm

A close-up of the Shallow Handpiece, showing a circular treatment area with a diameter of 5.2cm. The handpiece is grey with a green ring around the treatment area.

Deep Handpiece

Focal Depth 12mm

- Treatment area 6.5*6.5cm
- Applicable Areas: Body slimming, Thick fat areas such as the waist, abdomen, thighs, buttocks, and sides, as well as stubborn fat areas.

6.5cm

A close-up of the Deep Handpiece, showing a circular treatment area with a diameter of 6.5cm. The handpiece is grey with a green ring around the treatment area.

Laser Therapy

685nm Wavelength

- Treatment area 6.7*6.7cm
- 7pc laser lamp
- Applicable Areas: conditions such as frozen shoulder, cervical spondylosis, lumbar muscle strain, acupuncture point irradiation and Shallow fat removal.

6.7cm

A close-up of the Laser Therapy Handpiece, showing a circular treatment area with a diameter of 6.7cm. The handpiece is grey with a red ring around the treatment area, which contains seven red laser lamps.

Shallow & Deep Reduced Fat

2450 MHz
microwaves



Facial Firming



Fat Reduction



Cellulite



Body Firming



Shallow



Deep

Pocket Facial Firming



Focal Depth **3mm**



Full face lifting



Neck skin firming



Increase skin elasticity

Laser Physical Therapy

685nm red light , relieve pain



Inflammation and wound repair

- Effective treatment for inflammation, neurological and traumatic pain.



Physical Therapy

- For conditions such as frozen shoulder, cervical spondylosis, lumbar muscle strain and acupuncture point irradiation



Reduce fat

- Helps break down fat
- Skin beauty



Four handle sizes option

01

POCKET HANDPIECE

Treatment area 4.2*4.2cm



Facial
Firming



Improve
sagging



Facial
cellulite



02

SHALLOW HANDPIECE

Treatment area 5.2*5.2cm




Reduce arm fat



Double chin



Skin laxity



03

DEEP HANDPIECE

Treatment area 6.5*6.5cm




Remove
waist fat



Remove
back fat



Remove
thigh fat



04

LASER THERAPY HANDPIECE

Treatment area 6.7*6.7cm

685nm red light Therapy



Physical
Therapy



Acupuncture
point irradiation





Helps break
down fat




FREEZEWAVE

Reduced Fat & Skin Tightening & Laser Therapy




「2.45GHZ」 MICROWAVE FREEZEWAVE

80% of Freezewave energy reaches fat layer to destroy fat cells, 20% superficial heat retained in the dermis




Pocket Handpiece

Facial Firming




Shallow Handpiece

Cellulite | Skin Tightening



Deep Handpiece

Cellulite | Fat Reduction



Laser Therapy

Physical Therapy



4 in 1

Freezwave Function



Facial firming



Remove double chin



Remove arm fat



Remove abdomen fat



Remove back fat



Skin tightening



Frozen shoulder



lumbar muscle strain



Advantage

Microwaves at 2450 Mhz ✓
Frequency, 1000W high-power
output energy

Output Power : 12~82J ✓
(body); 6~40J (face)

12.1-inch touch screen, ✓
Intelligent system, easy to operate

Four working handles (Pocket, ✓
Shallow, Deep, Laser Therapy)

80% energy acts on fat, ✓
20% acts on surface heat effect

Function: Facial firming, dissolves fat, ✓
Skin tightening, Laser therapy:
frozen shoulder, lumbar strain,
fat breakdown



Pocket handpiece



Shallow handpiece



Deep handpiece



Laser Therapy

Instrument Use Case





Product Function



Facial firming



Remove arm fat



Remove back fat



Frozen shoulder



Remove double chin



Remove abdomen fat



Skin tightening



lumbar muscle strain

Product Advantage

- Microwaves at 2450 Mhz Frequency, 1000W high-power output energy
- Output Power : 12~82J (body); 6~40J (face)
- 12.1-inch touch screen, Intelligent system, easy to operate
- Four working handles (Pocket, Shallow, Deep, Laser Therapy)
- 80% energy acts on fat, 20% acts on surface heat effect
- Function: Facial firming, dissolves fat, Skin tightening, Laser therapy: frozen shoulder, lumbar strain, fat breakdown.

Before and After



China Care Medical Equipment Co., Ltd
Address: Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China web: www.chinacaremedical.com

4 in 1



Body Contouring & Laser Therapy Machine



Pocket
handpiece



Shallow
handpiece



Deep
handpiece



685nm Laser
Physical Therapy

